

# Portia by Bray vaginal ring pessaries

Vaginal ring pessaries for treatment of pelvic organ prolapse



## Category

Women's Health

## Product overview

- Pessaries are used to treat symptomatic uterine prolapse in women
- Discreet and non-surgical the pessaries are uncomplicated healthcare solutions available in a range of sizes.
- Flexible PVC available in 50-110mm diameter, 12.5mm wall thickness
- Rigid Polythene available in 50-100mm diameter, 7.5mm wall thickness
- Pessary sizing, fitting and replacement is carried out under medical guidance
- Latex free, does not contain metal and designed for comfort and ease of use
- CE marked medical device, compliant with Directive 93/42/EEC
- Quality Management System ISO 13485
- Manufactured and packed in a UK GMP controlled environment
- Regulated by MHRA
- Please contact our Customer Service team for all options and information



Ref: 366  
PVC 12.5mm



Ref: 365  
Polythene 7.5mm

Product Ref	Product description
365	Rigid Polythene, 50mm - 100mm 7.5mm wall thickness
366	Soft Flexible PVC (Vinyl), 50mm - 110mm 12.5mm wall thickness
P.5024	Vaginal Ring Pessary Sizing Guide



Ref: P.5024



Bray Group Ltd  
+44 (0)1367 240736  
sales@bray-healthcare.com  
www.bray-healthcare.com

